

Post-Op Care


3M™ Tegaderm™ + Pad Film Dressing with Non-Adherent Pad

A waterproof, absorbent dressing.



Suggested Applications:

- Post-surgical dressing - protects closed, clean surgical wounds and donor sites.
- Lacerations, abrasions and burns, small incisions and excisions
- Light to moderately draining wounds




Benefits:

- All-in-one sterile dressing - reduces application time
- A bacterial and viral* barrier - helps prevent external bacteria and contaminants from entering the site
- A waterproof barrier - patient may bathe or shower without removing dressing
- Conforms to body contours, providing comfort
- Unique, non-adherent absorbent pad
- Wear time up to 96 hours

* In vitro testing shows that the transparent film of Tegaderm™, Tegaderm HPT™ and Tegaderm™ CHG dressings provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.

Catalog No.	Dressing Size	Pad Size	Dressings/Box	Boxes/Case
3M™ Tegaderm™ + Pad Film Dressing with Non-Adherent Pad				
3582	2 in. x 2-3/4 in. 5 cm x 7 cm	1 in. x 1-1/2 in. 2,5 cm x 4 cm	50	4
3584	2-3/8 in. x 4 in. 6 cm x 10 cm	1 in. x 2-3/8 in. 2,5 cm x 6 cm	50	4
3586	3-1/2 in. x 4 in. 9 cm x 10 cm	1-3/4 in. x 2-3/8 in. 4,5 cm x 6 cm	25	4
3587 Oval	3-1/2 in. x 4-1/8 in. 9 cm x 10.5 cm	1-3/4 in. x 2-3/8 in. 4,5 cm x 6 cm	25	4
3588	6 in. x 6 in. 15 cm x 15 cm	4 in. x 4 in. 10 cm x 10 cm	25	4
3589	3-1/2 in. x 6 in. 9 cm x 15 cm	1-3/4 in. x 4 in. 4,5 cm x 10 cm	25	4
3590	3-1/2 in. x 8 in. 9 cm x 20 cm	1-3/4 in. x 6 in. 4,5 cm x 15 cm	25	4
3591	3-1/2 in. x 10 in. 9 cm x 25 cm	1-3/4 in. x 8 in. 4,5 cm x 20 cm	25	4
3593	3-1/2 in. x 13-3/4 in. 9 cm x 35 cm	1-3/4 in. x 11-3/4 in. 4,5 cm x 30 cm	25	4



Post-Op Care

3M™ Tegaderm™ + Pad Film Dressing with Non-Adherent Pad

Directions for Use:

Dressing Selection: Choose a dressing size large enough to provide a margin that adheres to dry healthy skin around the wound site.

Site preparation:

1. Prepare the site according to institution protocol.
2. Clipping of hair at the site may improve adhesion. Shaving is not recommended.
3. Skin should be clean and dry prior to application. Apply 3M™ Cavilon™ No Sting Barrier Film to protect peri-wound skin.
4. Allow all prep liquids to completely dry before applying the dressing to prevent skin irritation and ensure good adhesion.

Application:

1. Open the package and remove sterile dressing.

2.  Peel the paper liner from the paper-framed dressing, exposing the adhesive surface.

3.  Position the framed window over the wound site and apply dressing. Do not stretch the dressing during application as tension can cause skin trauma.

4.  Press the dressing into place. Remove the paper frame from the dressing while smoothing down the edges.

5.  Seal securely around wound site. Firmly smooth adhesive border to the skin. For added securement, seal the edges of the dressing with remaining 3M™ Cavilon™ No Sting Barrier Film.

Removal:

1. Gently grasp an edge and slowly peel the dressing from the skin in the direction of hair growth. If necessary, press a small piece of 3M™ Transpore™ Medical Tape onto a corner to lift an edge of the dressing.



3M Canada Company
3M Canada Company
P.O. Box 5757
London, ON N6A 4T1
Canada
1 800 364-3577

www.3M.com/ca/healthcare

3M, Cavilon, Tegaderm and Transpore are trademarks of 3M. Used under license in Canada. Please recycle. Printed in Canada. © 2010, 3M. All rights reserved. 1006-02156E